




KONICA MINOLTA



Dynamic Digital Radiography Analysis Workstation

# DI-X1


**RETHINK WHAT'S POSSIBLE**

Giving Shape to Ideas



## X-ray in motion

It is time to shift from static images to dynamic images



---

Dynamic Digital Radiography available in general radiography environment

Frame rate: 15fps


---

Exposure time: Up to 20 seconds

# Dynamic digital radiography - Create a new value in new ways

The Dynamic Digital Radiography system is composed of the Dynamic Digital Radiography Analysis Workstation “DI-X1”, the cassette type “AeroDR3 HD2”, and a general-purpose radiography system\*1.

This breakthrough system generates dynamic images by continuously irradiating a pulsed X-ray and displaying a series of frame-by-frame images.



# In pursuit of easy observation of “motions”


## Visualizing “motions”

The image processing we developed to improve visibility for static images is also provided for the dynamic images.

BS-MODE  
(Bone Suppression)



FE-MODE  
(Frequency Enhancement)




LM-MODE (Option)



Summary image


LM-MODE : Dynamic




## Quantifying “motion”

The motion of structures such as diaphragms can be quantified and displayed graphically. By quantifying motions, the system can provide new perspective to the evaluation of symptoms and functions.

DM-MODE  
(Specific component tracking process)




TD-MODE (Tracheal diameter)




Lung Field Area




Lung Field Measurement



Measurement (Signal Value Change, Position Tracking)



Signal Value Change




Position Tracking

## Extraction of signal value changes associated with "motion" of lung tissues


By capturing the signal value changes associated with the behavior of lung tissues such as alveoli and pulmonary blood vessels, and showing them on the display, the system expresses the slightest motions of the living body. This is expected to contribute to the information collection for a certain level of functional evaluation, even in simple X-ray studies in the screening stage before using advanced diagnosis devices such as CT and MRI.

### PL-MODE (Reference frame rate calculation process)

Extracts and displays the image density changes in the lung field associated with breathing.




Concentration changes in the lung field during breathing cycle



### PH2-MODE (Option)

Extracts and displays the signal value change in the lung field associated with the breathing of blood vessels




### PH-MODE (Cross-correlation calculation process)

Extracts and displays concentration changes in the lung field associated with the beating of blood vessels.



Image density changes in the lung field during cardiac cycle




### PH2-MODE (Option)



### Signal value ratio per lung Region

Calculate the signal value ratio of 6 areas;

$$\text{Signal value ratio} = \frac{\text{Sum of signal values in the area}}{\text{Sum of signal values for the entire lung field}}$$






Signal value ratio per lung Region-Summary

# In pursuit of easy observation of “motions”

## Measurement Summary(option)

In addition to multipurpose measurement results, excursion, lung field area, tracheal diameter, and measurement results can be summarized for viewing and DICOM PACS export.

6



— The average value is the measured average value of the control group indicated in the reference literature.  
— Measured value

# Specifications

## ■ Main Specifications of Dynamic Digital Radiography Analysis Workstation DI-X1

Product sales name: Image Diagnosis Workstation Konica Minolta DI-X1

	DI-X1 server
CPU	Intel Xeon® Processor E5-1650 v4 or higher
Memory	At least 64GB as standard
Operating system	Microsoft Windows10 IoT Enterprise 64bit
SSD or HDD	SSD: 512GB or more, and HDD: 4TB or more
PCI slot	At least 2 slots available (PCI-Express x16)
Ethernet	Gb Ethernet
USB	4 ports or more
Keyboard	Full-size keyboard (USB or PS/2)
Mouse	Optical wheel mouse (USB)
Monitor	23-inch TFT monitor, etc. with full HD (1920 x 1080) resolution or higher

	DI-X1 client
CPU	Intel Core i7 processor 3.40 GHz or higher
Memory	16GB or more
Operating system	Microsoft Windows10 IoT Enterprise 64bit
Storage	HDD: 2TB or more
PCI slot	At least 2 slots available (PCI-Express x 16)
Ethernet	Gb Ethernet
USB	4 ports or more
Keyboard	Full-size keyboard (USB or PS/2)
Mouse	Optical wheel mouse (USB)
Monitor	23-inch TFT monitor, etc. with full HD (1920 x 1080) resolution or higher

## ■ Main Specifications of AeroDR3 HD2

Product sales name: Digital Radiography SKR 3000

	AeroDR3 HD2 1417size (for dynamic imaging)
Scintillator	CsI
Detector size	14" x 17"
Dimension	384(W) x 460(D) x 15 (H) mm
Weight	2.6kg
Pixel size*2	400µm
Available grids	34 lines / cm, 40 lines / cm
DQE	56% (1mR, 1cycle / mm) 72% (1mR, 0cycle / mm)
Antibacterial	An inorganic antibacterial agent kneaded into the exterior material (complied JIS Z 2801 / ISO22196 test)
Communication*2	Dedicated wired Ethernet connection / Wireless LAN (IEEE802.11a / 802.11n compliant)

	AeroDR3 HD2 1717size (for dynamic imaging)
Scintillator	CsI
Detector size	17" x 17"
Dimension	460(W) x 460(D) x 15(H) mm
Weight	3.2kg
Pixel size*2	400µm
Available grids	34 lines / cm, 40 lines / cm
DQE	56% (1mR, 1cycle / mm) 72% (1mR, 0cycle / mm)
Antibacterial	An inorganic antibacterial agent kneaded into the exterior material (complied JIS Z 2801 / ISO22196 test)
Communication*2	Dedicated wired Ethernet connection / Wireless LAN (IEEE802.11a / 802.11n compliant)



# Specifications

## Console CS-7

Product sales name: Image Diagnosis Workstation CS-7

	CS-7H
Monitor	17-inch high viewing angle LCD multi-touch monitor
Image processing	Automatic grayscale processing Frequency processing (F processing) Equalization processing (E processing) Hybrid processing (HF processing/HE processing) Hybrid smoothing processing (HS processing) REALISM processing (RF processing/RE processing) REALISM smoothing processing (RS processing) Grid suppression processing Automatic irradiation field recognition processing
Image storage capacity	80GB or more
I/O port	Ethernet 10base-T/100base-TX/1000base-T USB 6 ports or more
DICOM support	DICOM Storage SCU DICOM Basic Grayscale Print Management SCU DICOM Modality Worklist Management SCU DICOM Modality Performed Procedure Step SCU
Image output gradation	12-bit (4,096 gradations)
Image output	Host: up to 4 ch. Printer: up to 2 ch.
Dimension/Weight	360(W) × 435(D) × 175(H) / 12 kg
Rated voltage	AC 100V-240V (single phase 50/60Hz)
Power consumption	109W or less
Main accessories	Control unit power cable, Operation unit power cable, Control unit - Operation unit connection cable, Keyboard x 1 / Mouse x 1
Main options	Magnetic card reader for patient ID registration, bar code reader, additional LAN card, sub-monitor adapter

## System Configuration Examples



★ Specifications are subject to change without notice to improve performance. ★ Hardware products are subject to change without notice to improve performance. ★ "Dynamic Digital Radiography Analysis Workstation DI-X1" and "DI-X1" are the names of "Image Diagnosis Workstation Konica Minolta DI-X1".  
★ "AeroDR3 HD2" is the name of "Digital Radiography SKR 3000". ★ The KONICA MINOLTA logo and symbol are registered trademarks of Konica Minolta, Inc. in Japan and/or other countries. ★ AeroDR is a registered trademark or trademark of Konica Minolta, Inc. in Japan and/or other countries.

